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Task	Description
1 CleanSyringe	PreClean the syringe with Standard Wash 1_MeOH
2 ITSPPrep_PAL3	Condition the ITSP cartridge with conditioning solvent 1_MeOH
3 ITSPPrep_PAL3	Condition the ITSP cartridge with conditioning solvent 2_H2O
4 ITSPPrep_PAL3	Load sample onto the ITSP cartridge
5 ITSPPrep_PAL3	Wash the ITSP cartridge with wash solvent 1_pH6.0 LSHCY
6 ITSPPrep_PAL3	Wash the ITSP cartridge with wash solvent 2_H2O
7 ITSPBlowOut_PAL3	Blow Out the ITSP cartridge
8 CleanSyringe	Clean the syringe with Wash1_MeOH
9 CleanSyringe	Clean the syringe with Wash1_MeCl2
10 ITSPElute_PAL3	Elute the sample from the ITSP cartridge to the eluate tray
11 ITSPBlowOut_PAL3	Blow Out the ITSP cartridge
12 CleanSyringe	Clean the syringe with Wash1_MeOH
13 WaitForStartSignal	Wait for GC to become ready
14 CleanSyringe	PreClean the syringe with Wash1_MeOH
15 Repeat	*Optional syringe sample clean
16 Transfer	Syringe sample clean
17 RepeatEnd	*Optional syringe sample clean
18 Injection	Get sample and inject
19 CleanSyringe	Post Clean the syringe with Wash1_MeOH
20 MoveToHome	

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